



Single- nucleus RNA Sequencing

Advancements in single cell capturing and individual library construction techniques, combined with high-throughput sequencing, have revolutionized gene expression studies on a cell-by-cell basis. This approach allows for a comprehensive analysis of complex cell populations, avoiding the masking of heterogeneity that occurs when taking the average of all cells. In cases where certain cells cannot be made into single-cell suspensions, nucleus extraction from tissues and prepared into single-nucleus suspension for single-cell sequencing is necessary. BMK offers the 10× Genomics Chromium™ based single-cell RNA sequencing service, which is widely utilized in disease-related research, like immune cell differentiation, tumor heterogeneity, and tissue development.

Technical Workflow



Service Workflow



Bioinformation list

➤ Data Quality Control

- Raw data quality score
- Statistics on sequencing output
- Library quality control: Data saturation, Gene expression quantification
- Statistics on single nuclei counts and read counts per nuclei

➤ Inner / Inter sample analysis

- Gene expression matrix based cell clustering
- **Gene expression analysis:** Gene expression quantification and distribution
- **Differential expression analysis:** marker gene identification, known cell type identification

- Protein interaction on marker genes of each cluster
- Hierarchical clustering on DEGs
- Prediction on TF of DEGs
- Protein interaction analysis of DEGs
- Annotation of DEGs on disease database
- Function annotation and enrichment of DEGs and gene set of each cluster

Service Advantages

Single-nucleus vs Single-cell

Tissue not suitable for single cell suspension preparation

Single-nucleus	Single-cell
Unlimited cell diameter	Cell diameter: 10-40 μ m
The material can be frozen tissue	The material must be fresh tissue
Low stress of frozen cells	Enzyme treatment may cause cell stress reaction
No red blood cells need to be removed	Red blood cells need to be removed
Nuclear expresses bioinformation	The whole cell expresses bioinformation

Cell / Tissue	Reason
Unfresh frozen tissue	Unable to get fresh or long-saved organizations
Muscle cell, Megakaryocyte, Fat...	Cell diameter is too large to enter the instrument
Liver...	Too fragile to break, unable to distinguish single cells
Neuron cell, Brain...	More sensitive, easy to stress, will change the sequencing results
Pancreas, Thyroid...	Rich in endogenous enzymes, affecting the production of single cell suspension



Service Specifications

Sample Requirements


Library	Platform	Data volume
10 \times Genomics single-nucleus library	10 \times Genomics Illumina PE150	100,000 reads/cell approx. 100-200 Gb

Cell	Tissue
Cell number: $>2 \times 10^5$ Cell conc. at 700-1,200 cell/ μ L	≥ 200 mg


Demo Results




Heatmap: gene expression clustering across cell sub-populations




Cell trajectory analysis/pseudotime



Cell sub-population identification



Cell cycle identification of sub-population



Featured Publications

Year	Journal	Title
2023	Int J Biol Sci	Integrating Spatial Transcriptomics and Single-nucleus RNA Sequencing Reveals the Potential Therapeutic Strategies for Uterine Leiomyoma
2023	Phytomedicine	Qi-Po-Sheng-Mai granule ameliorates Ach-CaCl ₂ -induced atrial fibrillation by regulating calcium homeostasis in cardiomyocytes



Biomarker Technologies (BMKGENE) GmbH

BioZ, Johann-Krane Weg tech@bmklcloud.com
 42, 48149 Münster, Germany www.bmkgene.com

Copyright©2009-2023 Biomarker Technologies (BMK) GmbH.
 All Rights Reserved. Information and specifications are subject
 to change at any time without notice.