


# Whole Transcriptome Sequencing


## Product Introduction

Whole transcriptome sequencing is designed to profile all types of RNA molecules, including coding (mRNA) and non-coding RNAs (including lncRNA, circRNA, and miRNA) which are transcribed by specific cells at a certain time. Taking advantage of NGS technology, whole transcriptome sequencing aims to reveal comprehensive regulatory networks of circRNA-miRNA-mRNA based ceRNA and perform joint RNA analysis.

## Technical Features



## Service Workflow



## Bioinformation List

- ▶ RNA expression overview;
- ▶ Differential expression overview;
- ▶ Co-expression analysis;
- ▶ Competing endogenous RNA(ceRNA) network;
- ▶ Integrated pathway analysis on key genes;
- ▶ Targeting relations between differential RNAs.

## Service Advantages

- ▶ Estimation of all types of RNAs in terms of counts, expression, and chromosome-based relative expression.
- ▶ Identification of differentially expressed RNAs and corresponding expression.
- ▶ Gene co-expression analysis and ceRNA network analysis.
- ▶ BMKGENE has completed 2000+ whole transcriptome sequencing projects.
- ▶ After-sale services: After-sale services are valid for 3 months upon project completion, including project follow-up, trouble-shooting, results Q&A, etc.

## Service Specifications


Library	Platform	Recommended data	Data QC	Time
Total RNA	Illumina PE150&SE50	Circ/lnc/mRNA: 16G; miRNA: 10M reads	Q30≥85%	40 days

## Sample Requirements


Amount	Qubit Conc. (ng/μL)	Integrity
Conc. ≥100 ng/μL; Volume ≥ 10 μL; Total ≥ 1 μg	OD260/280=1.7-2.5 OD260/230=0.5-2.5 Limited or no protein or DNA contamination shown on gel.	For plants: RIN≥6.5; For animals: RIN≥7.0; 5.0≥28S/18S≥1.0; limited or no baseline elevation

## Demo-Results

1. Estimation on all types of RNAs based relative expression



2. Integrated KEGG pathway network



3. ceRNA network analysis



## Featured Publications

Year	Journal	Paper	Article	DOI
2023	Pharmacological Research	Targeting hnRNPC suppresses thyroid follicular epithelial cell apoptosis and necroptosis through m6A-modified ATF4 in autoimmune thyroid disease	Disease treatment	10.1016/j.phrs.2023.106933
2023	Aging	Whole-transcriptome sequencing analysis reveal mechanisms of Yiqi Huoxue Yangyin (YHY) decoction in ameliorating D-gal-induced cardiac aging	Disease treatment	10.18632/aging.204532
2022	Frontiers in Plant Science	Identification and validation of coding and non-coding RNAs involved in high-temperature-mediated seed dormancy in common wheat	Environmental adaptation	10.3389/fpls.2023.1107277



### Biomarker Technologies (BMKGENE) GmbH

BioZ, Johann-Krane Weg 42, 48149 Münster, Germany
 tech@bmkgene.com
 www.bmkgene.com

Copyright©2009-2023 Biomarker Technologies (BMK) GmbH.  
All Rights Reserved. Information and specifications are subject to change at any time without notice.